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Quality and Reliability

Simplex Rapid Powder and Liquid

Made in England by:
Associated Dental Products Ltd
Purton, Swindon,
Wiltshire, SN5 4HT UK
Tel: +44 (0)1793 770256
e-mail: sales@kemdent.co.uk
www.kemdent.co.uk

INSTRUCTIONS FOR USE (ENGLISH)

INDICATIONS

Mix Simplex Rapid Powder and Liquid to form a selfcure (autopolymerising) acrylic resin. Simplex Rapid conforms to BS EN ISO 20795-1 Type II Class I for repairs, indirect relines and additions to removable prosthetic appliances. Will bond to acrylic teeth when cured by Hydroflask.

PROPORTIONS AND MIXING

Dispense 10 drops of liquid into a mixing vessel and sprinkle powder onto the liquid whilst vibrating the dish until all the liquid is absorbed. Shake off any excess powder. Add 5 drops of liquid and stir slowly with a metal spatula for 10 seconds (45 seconds if room temperature is below 18°C) The mix should now be fluid, giving continuous flow when poured. Apply immediately.

REPAIRS

Cast the plaster model as usual. Chamfer the broken edges of the denture. If necessary protect both the fitting surface of the denture and the teeth adjacent to the repair with petroleum jelly. Coat model with plaster separating solution (Kemdent Liquid Tin Foil). Seat the denture on the model and secure with Kemdent Sticky Wax if necessary. Flow the fluid mix into the fracture. (Do not use after mix ceases to flow freely).

LABORATORY RELINING

Cast the plaster model as usual. Produce locating matrix for denture in the usual way. Reduce fitting the surface of denture. Coat model with plaster separating solution (Kemdent Liquid Tin Foil). Make an appropriate size mix of Simplex. Spread onto relined area of denture. When surface dries locate on model using matrix. Seat with firm pressure for 10 seconds. Secure firmly and place in Hydroflask.

Other applications

Replacing and adding teeth
Addition of clasps

CURING

For best results cure in hydroflask at 40°C for 10-20 minutes at 30psi
Remove carefully, trim and polish

Will also self cure in 20 minutes

NOTE: Prior to curing, ensure model has not been allowed to dry out as trapped air may affect the quality of acrylic.

NEGATIVE INTERACTIONS WITH DENTAL MATERIALS: None noted

CONTRAINDICATIONS AND SIDE EFFECTS

The curing of Simplex involves an exothermic reaction so care must be taken if it is handled whilst setting. Cured material may cause an allergic reaction in those susceptible Methyl Methacrylate. May cause sensitisation by skin contact.

PRECAUTIONS

For use by the Dental Professional only

Powder: replace cap after use
Do not inhale dust
Store between 5-25°C

Liquid: Mixture based Methyl Methacrylate Monomer

RISK PHRASES

Highly flammable Flash Point below 11°C
Harmful by inhalation, in contact with skin and if swallowed.
Danger of cumulative effects
Irritating to eyes, respiratory system and skin
May cause sensitisation by skin contact

SAFETY PHRASES

Keep container in well ventilated place
Keep away from sources of ignition - no smoking
Do not empty in to drains
Take precautionary measures against static discharge

Avoid prolonged exposure to sunlight

In case of contact with eyes and skin: flush with copious volumes of water.
Seek medical attention
Store below 25°C. Replace cap immediately after use.

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