

ASSOCIATED DENTAL PRODUCTS LTD
MATERIAL HEALTH AND SAFETY DATA SHEET

Version: 7

Updated on: 06/10/05

1. IDENTIFICATION	
1.1: Product Name And Code: LAB001 Zonalin Powder 1.2: Company Identification: Associated Dental Products Ltd Kemdent Works Purton Swindon Wiltshire SN5 4HT Tel: +44 (0) 1793 770 256 Fax: +44 (0) 1793 772 256 1.3: Emergency Tel: +44 (0) 1793 770 256	
2. COMPOSITION	3. HAZARDS IDENTIFICATION
2.1: Zinc oxide BP & zinc acetate GBR 2.2: Hazardous Ingredients: 2.3: 2.4: R Phrases: For full text of R Phrases see section 16	None classified
4. FIRST AID MEASURES	5. FIRE FIGHTING MEASURES
Eyes: Flush with plenty of water Skin: Wash and flush with plenty of water Mouth: Rinse mouth with water Inhalation: Remove to fresh air	Extinguishing Media: Non flammable Exposure Risk: Normal protective clothing
6. ACCIDENTAL RELEASE MEASURE	7. HANDLING AND STORAGE
Precautions: None specified	7.1: Handling: No special precautions 7.2: Storage: Keep container closed 7.3: Specific Uses:
8. EXPOSURE CONTROLS/PERSONAL PROTECTION	9. PHYSICAL AND CHEMICAL PROPERTIES
8.1: Exposure Limit Values: 8.2: Exposure Controls: The following information is given as general guidance Respiratory Protection: Use in a well ventilated area Hand Protection: Use gloves if considered necessary Eye Protection: None Skin Protection: Normal overalls are adequate	Appearance: Slightly creamy free flowing powder Odour: None PH: Melting Point: 1970C Boiling Point: Flash Point: Explosive Properties: Oxidizing Properties: Vapour Pressure: Relative Density: 5.6 SG Solubility In Water: Insoluble Solubility In Organic Solvent: NA Partition Coefficient: Other Data:
10. STABILITY AND REACTIVITY	11. TOXICOLOGICAL INFORMATION
Hazardous Reactions: 10.1: Conditions To Avoid: Release of dust 10.2: Materials To Avoid: None 10.3: Hazardous Decomposition Products:	Inhalation: Non Toxic Skin Contact: Eye Contact: Ingestion: Long Term Exposure:
12. ECOLOGICAL INFORMATION	13. DISPOSAL CONSIDERATIONS
12.1: Ecotoxicity: 12.2: Mobility: 12.3: Persistence and Degradability: 12.4: Bioaccumulative: 12.5: Other Adverse Effects:	Dispose of according to local regulations
14. TRANSPORT INFORMATION	15. REGULATORY INFORMATION
No special precautions	Hazard Symbol: Safety Phrases:
16. OTHER INFORMATION	
16.1: Shelf Life: 2 Years	