



Introducing PlasterSafe – Plaster Solvent & Tray Cleaner

Kemdent is celebrating the release of PlasterSafe, our new plaster solvent and tray cleaner. PlasterSafe is alcohol free, ready to use solution that aids the safe removal of dental plaster and stone from the fitted surface of dentures and the hard to reach areas.

PlasterSafe is kind to technician's hands, it is a water based solution with no harmful chemicals. It is also easy to use, just remove any excess plaster or alginate, then immerse your trays in PlasterSafe. This can either be done rapidly in an Ultrasonic bath or soaked for as long as possible. The plaster dissolves and the trays can be rinsed quickly and efficiently in clean water.

PlasterSafe is cost effective, saves time and effort, it is the ideal safe solution for the busy dental laboratory.



For further information on PlasterSafe contact Jenny on +44 (0)1793 770256
E-mail jenny@kemdent.co.uk or visit our website www.kemdent.co.uk